



# HOW TO ENSURE YOU GET THE FULL BENEFIT FROM YOUR JUSTFLŌ NEBULISER.

justflō  
flōcanix®




## Prepare the device and font

Take the refill font out from the packaging by tearing the foil. Ensure the battery device is clean and free of moisture.



## Attach the font to the device

Slide in the font. You will notice the font attaches itself to the device, which means the font and device are connected through a magnetic lock.



## Inhale and exhale slowly

Draw on the mouthpiece until your mouth fills with the aerosol, then inhale and exhale slowly. No buttons needed.